East West Immunogenetics Conference Prague, Czech Republic March 5 – 6, 2009

## Practical Aspects and Typing Strategy in Molecular ABO and RH Blood Group Diagnosis

Martina Prager prager.martina@bag-healthcare.com





- General and specific application
- Clinical cases
- Investigation strategy for clarification of unclear serological blood group typing
- Examples of genomic nomenclature (ABO and RHD)
- Description of the BAGene DNA SSP Kits



# **General Application**

- Discrepant, weak or unexpected results in blood group serology
- Donor Pass
- A Rhesusprophylaxis / Transfusions
- Today a donor, tomorrow a patient and vice versa



# **Specific Application**

- Genotype multi-transfused recipients
- Genotype patients after ABO-incompatible bone marrow transplantation
- Determine RHD zygosity of partners of alloimmunized D negative women before pregnancies
- Genotype Rh D negative donors with C or E in order to exclude the presence of the RHD gene and thus preventing anti-D alloimmunization of recipients caused by hidden Rh D variants in RBC units
- Identify genotype in case of weakly expressed Rh D (e.g. DEL) in donors
- Confirm weak D genotype in recipients in order to avoid the donation of Rh D negative blood units
- Quality control of serological methods
- External Quality Assurance trials



# Molecular Genetic Blood Group Typing

# Principle of PolymeraseChainReaction using SequenceSpecificPrimers

Perfect match→ Amplification (specific Allel) Mismatch → no Amplification (unspecific Allel)

**Maniatis et al.** 1989. Molecular Cloning: A Laboratory Manual. New York: Cold Spring Harbour Laboratory

Beutler E et al. 1990. BioTechniques 9:166

Olerup O, Zetterquist H 1992. Tissue Antigens 39:225-235





### **Investigation Strategy** ABO blood group typing **Serological Determination** unambiguous determination ABO Weak antigen expression or of ABO including subgroups unclear weak and unexpected reactions in reverse testing respectively **BAGene BAGene ABO-TYPE variant ABO-TYPE** confirmation by ABO genotype molecular genetics clarification by sequencing in case of ambiguous genotype

**BAG** HEALTH CARE

# **Investigation Strategy**


### A Rhesus typing






Result according to Transfusion guideline

Rhesus D positive CcD.ee Kell positive










microwell plate antisera dilution 1:5


agglutination strength (1) -1

tube test immediate spin

questionable agglutination with monoclonal Anti-D


**BAG** HEALTH CARE







Worksheet und Auswertetabelle / Worksheet and Evaluation diagram




# Final result weak D type 1

# Rhesus prophylaxis not required


Transfusion of D positive red blood cells is possible






Result according to Transfusion guideline

Rhesus D negative, Ccddee Kell negative









# microwell plate antisera dilution 1:5

negative reaction





### **BAG** HEALTH CARE


BAGene	CE	LOT 0902 WD
	IVD	2010-08
Weak D-TYPE	Ĩ	REF 6647

Worksheet und Auswertetabelle / Worksheet and Evaluation diagram

	Reaktions-Nr. / Reacti	on No.	1	2	3		4	5	6	7	8
	PCR-Produkt (Größe i PCR product (size in bp)	150	126	16	5	101	130 83	112	198 83	153	
				wea	k D Al	lele /	weak D a	lleles			
	weak D type 1		+	-	-		-	-	-	-	-
	weak D type 2		-	+	-		-	-	-	-	-
$\Rightarrow$	weak D type 3		-	-	+		-	-	-	-	-
	weak D type 4.0, 4.1		-	-	-		+	-	-	-	-
	weak D type 4.2, DAR	-	-	-		+	130	-	-	-	
	weak D type 5	-	-	-		-	-	+	-	-	
	weak D type 11 (haploty	oe cDe)	-	-	-	,	-	-	-	198	-
	RHD(M295I) (haplotype (	CD <sub>el</sub> e)	-	-	-		-	-	-	198	-
⇒	weak D type 15		-	-	-		-	-	-	-	+
	weak D type 17		-	-	-		-	83	-	83	-
	weak D type 4.2, 17		-	-	-		+	130 83	-	83	-
	Weak D type 11 / RHD(M2	95I), 17	-	-	-		-	83	-	198 83	-
	RHD pos. oder / or RHD r	-	-	-		-	-	-	-	-	
				1							
	Ergebnis Result	1	2	3	4	5	6	7	8	Geno Genot	

Ergebnis 1 Result	-	5	7	5	0	'	0	Genotyp

Proben-ID / Sample-ID: 0627 9624 Name: Geb.-Datum / Birthdate: Ergebnis/Result: => D. partial Datum/Date: 25. 10.04 Unterschrift / Signature: Lesames



### **BAGene**

**Partial D-TYPE** 




**BAG**HEALTH CARE

#### Worksheet und Auswertetabelle / Worksheet and Evaluation diagram

Reaktions-Nr.						1	2	3	4	5	6	7	8	9	10	11	12	13	14	1
PCR-Produkt (Größe in bp) PCR product (size in bp)						134	146	118	135	132	132	120	673	184	132	166	107	117	129	14
RHD Exons * Mix 8: Intron 7/Ex	type	D,		D <sub>3</sub>	D4	Ds	D <sub>6</sub>	D,	D <sub>7/8</sub> *	D <sub>9</sub>	D <sub>10</sub>	D <sub>2</sub>	D <sub>6</sub>	D,	D,	C				
Spezifität / Specif					.,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,		-	-	-	-		S	tanda	rd RI	HD			-	_	-
D-positive			+			+	+	+	+	+	+	+	+	+	+	-	-	-	-	Ι.
b positive			+			-							HD ne							-
D-negative, cc			+			-	-	- 1	-	-	-	-		-gati	-	-	-	-	-	Ϊ.
D-negative, Cc oder / or	cc		-		-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	t.
D negutite, ce ouer, or			+			-		8			RHD	Vari	anten	/ RH	D Var	riants				-
D cat, II at				CD	e	+	+	+	+	+	+	+	+	+	+	-	-	-	+	
D cat.IIIa, IIIc, III type 4	4		-	De, CDe	e, CDe	+	+	-	+	+	+	+	+	+	+	-	-	-	-	
D cat, IIIb **				cDe	2	+	-	+	+	+	+	+	+	+	+	-	-	-	-	•
D cat. IVa				cDe	3	+	+	-	+	+	+	-	-	+	+	-	-	-	-	1.
D cat. IVb, IVb(J) nt				CD	e	+	+	+	+	+	+	-	-	-	+	-	-	-	-	
D cat. IV type 3				CDe	e	+	+	+	+	+	-	-	-	-	+	-	-	-	-	
D cat. IV type 4				CD	e	+	+	+	+	+	+	-	+	+	+	-	-	-	+	•
D cat. Va, Va-like, Va-ass	ociate	d, DE	8S	CDe	e	+	+	+	+	-	+	+	+	+	+	-	-	-	-	
D cat. VI type 1						+	+	+	-	-	+	+	+	+	+	-	-	-	-	
Cat. VI type 2				CDe	2	+	+	+	-	-	-	+	+	+	+	-	-	-	-	
D cat. VI type 3				CDe	e	+	+	-	-	-	-	+	+	+	+	-	-	-	-	
D cat. VI type 4 [DHMii]	n.t.			CDe, C	DE	+	+		-	-	+	+	+	+	+	-	-	-	-	
D cat. VII				CDe	e	+	+	+	+	+	+	+	+	+	+	+	-	-	-	•
DAR (weak D-Type 4.2)				cDe		+	+	+	+	+	+	+	-	+	+	-	-	-	-	-
DAU				cDe		+	+	+	+	+	+	+	+	+	+	-	-	-	-	
DBT type 1						+	+	+	+	-	-	-	-	+	+	-	-	-	-	
DBT type 2 nt.				CDe	9	+	+	+	+	-	-	-	-	-	+	-	-	-	-	•
DFR + oder / or RHD psi						+	+	+		+	+	+	+	+	+	-	-	-	-	-
DHMi				cDE		+	+	+	+	+	+	+	+	+	+	-	+	-	-	
DNB				CDe	2	+	+	+	+	+	+	+	+	+	+	-	-	+	-	
RHCE-D(5)-CE (DHAR (R	133))			cDe		-	-	-	-	+	-	-	-	-	-	-	-	-	-	
RHCE(1-9)-D(10)				cdE		-	-	-	-	-	-	-	-	-	+	-	-	-	-	•
RHD-CE(2-9)-D			_	Cde		+	+	-	-	-	-	-	-	-	+	-	-	-	-	· ·
RHD-CE(3-7)-D d(C)c	e'			Cde		+	+	-	-	-	-	-	+	+	+	-	-	-	-	•
RHD-CE(4-7)-D			-	cdE		+	+	+	-	-	-	-	-	+	+	-	-	-	-	
RHD-CE(8-9)-D			_	Cde		+	+	+	+	+	+	+	-	-	+	-	-	-	-	•
RHD(delEx9)				CDe	_	+	+	+	+	+	+	+	+	-	+	-	-	-	-	-
RHD(Q41X) nt.				Cde		-	+	+	+	+	+	+	+	+	+	-	-	-	-	•
DCS, [DFW, DHR, DIM, DNU] <sup>n.t.</sup> diverse / various				+	+	+	+	+	+	+	+	+	+	-	-	-	-			
Ergebnis	1	-	-			area a		an 140	ted cu	a lana		12	14 4	- 1	Phän	otyp	e 1	Ge	noty	р
Result	1	2	3	4	5	6	7	8 9	) 10	11	12	13	14 1		Pheno				notyp	

0627 9624 Proben-ID / Sample-ID: Name: Geb.-Datum / Birthdate: Ergebnis / Result: D. VI Datum / Date: 26.10.04 Unterschrift / signature: Watue



									BAGen	е			HEALTH CARE		
E	BAGe <u>n</u> e	9				CE		LOT 0902 WD		Partial D-	ГҮРЕ	DFR (serologisc (hemi- bzw. ho Bei fehlender Ir mit BAGene R werden.	e Bande in Reaktion Nr. 4 weist h schwach D-positiv) oder RHD morzygot, serologisch D-negativ) formation zur Serologie kann RHD H-YYE bestätigt oder ausgeschloi nd in reaction no. 4 may indicate i tith anti-D) or RHD psi (hemi- or ho	psi CE hin. psi ssen IVD	LOT 0812 PD
	Veak D-	in	a		re	SU	llt	: D	cat	egory	VI	ty	/pe		diagram           13         14         15           117         129         140           D <sub>7</sub> D <sub>7</sub> D <sub>8</sub>
PCR-Produkt (Grö PCR product (size in weak D type 1 Rhesus prophylaxis required															
	weak D type 2 weak D type 3 weak D type 4.0, 4.1 weak D type 5 weak D type 11 (hap) weak D type 11 (hap)														
₩ ₩	RHD(M295I) (haploty veak D type 15 veak D type 17 veak D type 4.2, 17	0	0	d	Ce		S	is re	equ						
	Veak D type 11 / <i>RHD</i> (M295I), 1 RHD pos. oder / or <i>RHD</i> neg.	7 – –	-	-	-	83 -	-	198 83 – – –		RHCE-D(5)-CE (DHAR (Rh33)) RHCE(1-9)-D(10) RHD-CE(2-9)-D RHD-CE(3-7)-D d(C)ce <sup>5</sup>	CDe cDe cdE Cde Cde	+ + + + + +		· · · · · ·	
E	Ergebnis 1 Result	2	3	4	5 6	7	8	Genotyp Genotype	]	RHD-CE(4-7)-D RHD-CE(4-7)-D RHD(delEx3) RHD(Q41X) <sup>AL</sup> DCS, [DRW, DHR, DIM, DNU] <sup>AL</sup> Ergebnis 1		+ +	+ + + + + +	- + + - + - + - + - + + + + - + + + +	 
Proben-ID / Sam	nple-1D: 0627 968	24		٦٢		Ge	lbild / Ge	l Picture		Result				Phenotype	e Genotype
Name:					-	Dw	0627	9624	Probe	en-ID / Sample-ID: 0627 96	24			Gelbild / Gel	
GebDatum / B									Name					kart. 0627	9624
Ergebnis / Resul	t: => D. part	ial					-		-	Datum / Birthdate:					
Datum / Date:	25.10.04									mis/Result: D. <u>VI</u> m/Date: 26.10.04				-	
Unterschrift / Si	ignature: lesques							10 V 1		schrift/signature: Usatue					

## **Genomic Nomenclature - ABO blood groups**

### \* The Blood Group Antigen Gene Mutation Database

http://www.ncbi.nlm.nih.gov/projects/gv/mhc/xslcgi.cgi?cmd=bgmut/systems\_alleles&system=abo

ABO Phenotype	ABO Genotype					
	"old"	* BGMUT				
A <sub>1</sub>	A1	ABO*A101				
A <sub>2</sub>	A <sup>2</sup>	ABO*A201				
<b>A</b> <sub>3</sub>	A <sup>3</sup>	ABO*A301				
A <sub>el</sub>	A <sup>el</sup>	ABO*Ael01				
Aw	A <sup>w</sup>	ABO*Aw04, 06, 07, 11				
A <sub>x</sub>	A <sup>x</sup>	ABO*Ax01, 02, 03, 04, 05, 06				
В	B <sup>1</sup>	ABO*B101				
<b>B</b> 3	B <sup>3</sup>	ABO*B302				
Bw	$B^{w}$	ABO*Bw09				
Bx	B <sup>x</sup>	ABO*Bx01				
0	0 <sup>1</sup>	ABO*001				
0	0 <sup>1</sup> v	ABO*002				
0	0 <sup>2</sup>	ABO*003				



### **Genomic Nomenclature – Rhesus D**

### The Rhesus Site (RHD Mutation Database)

http://www.uni-ulm.de/%7Efwagner/RH/RB/

RhD Phenotype	RHD Genotype (Allele)	Mutation
D positive	RHD	./.
D negative	RHD deletion	Deletion of RHD
partial D, DVI	DVI	Hybrid
partial D, DFR	DFR	Hybrid
partial D, DBT	DBT	Hybrid
partial D, Rh33	DHAR	Hybrid
D negative, (C)cde <sup>s</sup>	RHD-CE(3-7)-D	Hybrid
partial D, DVII	DVII	Single nucleotide exchange
weak D	e.g. weak D type 1, 2, 3	Single nucleotide exchange
partial D	DHMi	Single nucleotide exchange
partial D, DNB	DNB	Single nucleotide exchange
weak D	weak D type 4.0/4.1, DAR/weak D type 4.2	Multiple nucleotide exchange
partial D, DAU	DAU	Multiple nucleotide exchange
D <sub>el</sub>	DEL types RHD(K409K), RHD(M295I)	Splice site
D negative	RHD(W16X)	Stop codon
D negative	RHDΨ	Multiple nucleotide exchange, deletion and insertion

**BAG** HEALTH CARE




### BAGene

SSP kits for determination of

- ABO blood groups
- RH types
- Kell, Kidd, Duffy systems
- MNS system
- HPA specificities

on a molecular genetic basis



# **Contents of our BAGene SSP kits**

PCR plates or strips with prealiquoted, dried and colored reaction mixes containing allele specific primers, internal control primers (specific for the HGH gene) and nucleotides.

10 x PCR buffer

PCR strip caps

Worksheets and Evaluation Diagrams

Instructions for Use



### **RHD-TYPE** Asia

REF 6649

CE

Kit zur molekulargenetischen RHD-Bestimmung einschließlich schwacher und selt RHD-Varianten (Weak D Typ 15, 17, RHD(K409K)) Kit for determination of RHD including weak and rare RHD variants (weak RHD(K409K)) on a molecular genetic sis

a Acone ONA-SSP Lites

and to a stand

a tago tot

sup sher S.Seasserge

Ser. spor and the state of the state

\* sonton

South Bran

Selen stor Service State Constanting of the and alight

Patra Pres - Patrice

- the Direct \*Disserver

P.O.T.

FORM

Hadre - ALANS

1. Contraction AN' PHAN Loss miles PCR-Streifon (M CR strips (Mix 1-7) gebrauchefet \* \* \* \* / ready to use, prealigu


iffer

caps ns for use 'iagrams

614

### **Evaluation of Results - RHD-RHCE genotyping**

### **BAGene**





 Amtsgerichtisstraße 1-5
 Tel.: +49 (0) 6404 / 925-0
 www.bag-healthcare.com

 35423 Lich/Germany
 Fax: +49 (0) 6404 / 925-250
 info@bag-healthcare.com

 Auftrags annahme/Ordering:
 Customer Service:

 Tel.: +49 (0) 6404 /925 -450
 Tel.: +49 (0) 6404 /925 -125

 Fax: +49 (0) 6404 /925 -460
 Fax: +49 (0) 6404 /925 -421

 verkau@bag-healthcare.com
 service@bag-healthcare.com

#### Kurzanleitung / Short Instructions

 Die gewünschte Anzahl BAGene-Platten/-Streifen aus dem Gefrierschrank (-20...-80°C) nehmen und den 10 x PCR-Puffer bei Raumtemperatur auftauen.

Remove the required number of BAGene-plates/-strips from -20...-80°C and thaw the 10 x PCR-buffer.

2. Der erste Reaktionsmix ist durch Aufdruck der Lot-Nr. markiert.

The first reaction mix is marked by the printed lot number.

 Den Mastermix, bestehend aus 10 x PCR-Puffer, DNA-Lösung, Taq-Polymerase und Aqua dest. zusammenpipettiren und gründlich vortexen. Die verschiedenen BAGene DNA-SSP Kits werden mit dem gleichen Mastermix angesetzt und sind daher miteinander kombinierbar.

Pipet the mastermix, consisting of 10 x PCR-buffer, DNA-solution, Taq-Polymerase and aqua dest. and mix well. The different BAGene DNA-SSP Kits do all work with the same mastermix and can therefore be combined.

#### Zusammensetzung des Mastermixes in Abhängigkeit von der Anzahl der Reaktionsmixe Composition of the Mastermix depending on the number of reaction mixes

No. of mixes	Aqua dest.	10 x PCR-buffer	DNA-sol. (50-100 ng/µl)	Taq-Polymerase (5 U/µl)	total volume app.	
1	8	1	1	0,08	10	μ
2	16	2	2	0,2	20	μ
6*	50	7	7	0,5	65	μ
7	70	9	9	0,7	90	H
8	80	10	10	0,8	100	H
9	88	11	11	0,9	110	H
10	96	12	12	1,0	120	F
11	104	13	13	1,0	130	ŀ
12	112	14	14	1,1	140	F
13	128	16	16	1,3	160	L
14	136	17	17	1,4	170	H
15	144	18	18	1,4	180	ŀ
16	152	19	19	1,5	190	H

- Bei abweichenden DNA-Konzentrationen sind die Mengen von DNA und Wasser entsprechend zu variieren. For different DNA concentrations, the quantities of DNA and water must be varied accordingly.
- Mastermix f
  ür 6 Reaktionsmixe wird wegen des geringen Volumens Taq-Polymerase als Mindestansatz empfohlen.

Minimum preparation of mastermix for 6 reaction mixes is recommended due to the small volume of Taq-Polymerase.

 Nach dem Vortexen werden von diesem Gemisch umgehend je 10 µl zu den angetrockneten Reaktionsmixen pipettiert.

After vortexing add 10 µl of this mixture immediately to the dried reaction mixes.

 Die Reaktionsgefäße werden mit den dafür vorgesehenen Deckeln dicht verschlossen. Die Platten/Streifen werden leicht bewegt, um das Pellet auf dem Gefäßboden etwas anzulösen. Die gesamte Reaktionslösung soll sich im Gefäßboden befinden.


Tightly close the tubes with the respective strip caps. Slightly move the plates/strips to dissolve the pellet at the bottom of the tube. All PCR-solution should settle on the bottom of the tube.


- Nach der PCR und Auftrennung der Amplifikate im Gel erfolgt die Auswertung. After PCR and separation of the amplicons in the gel, the evaluation can be carried out.
- 7. Alle BAGene Reagenzien sind nach Gebrauch wieder bei -20...-80°C zu lagern.

Return all BAGene reagents to -20...-80°C after use.

> Ablesen der Ergebnisse in Pfeilrichtung / Read the results in direction of the arrow


Beispiel: Gelbilder der Reaktionen 1 und 2







#### **ABO** Genotyping




BAGene ABO-TYPE variant genotype O<sup>1</sup> A<sup>x</sup>

#### **Genotyping partial D**



BAGene Partial D-TYPE genotype D cat. VI type II




BAGene D Zygosity-TYPE

#### Genotyping MNS



BAGene MNS-TYPE

### RHD, RHCE Genotyping



BAGene RH-TYPE genotype RHD(K409K) Ccee

### Genotyping weak D



BAGene Weak D-TYPE genotype weak D type 4.2

#### Genotyping KEL, JK, FY



BAGene KKD-TYPE genotype KEL\*1/KEL\*2; JK\*A/JK\*B; FY\*B/FY\*null01

#### Genotyping HPA





BAGene HPA-TYPE genotype HPA-1a/a; 2a/b; 3a/a; 4a/a; 5a/b; 15 a/b

### **Intended use**

The molecular determination of blood group antigens using SSP Kits has to be performed after serology

We offer these assays as a supplementary technique to investigate weak or discrepant serological findings. The current assays are not intended to replace serology.

In case of discrepant or unclear genotyping results, transfusion guidelines have to be followed in accordance with serological typings. Final clarification by sequencing analysis is recommended.



### Conclusion

SSP technique is helpful to resolve most of the problems caused by discrepant or doubtful serological results

**Easy to handle** 

Questionable cases in donor, recipient and patient typing can be examined with acceptable cost.

Not suitable for high throughput

Not suitable to replace serology



### References

 Frequency of DAU, RHD psi, Cde<sup>s</sup> in South African blood donors. Maas JH, Binder E, Smart E, Prager M, Legler T. Transfus Med Hemother 2004;31(suppl 3):PS412

RHD and RHCE Genotyping in South African Blood Donors with Prepipetted PCR-SSP Kits.

Legler TJ, Binder E, Smart E, Prager M, Maas JH. Transfusion 2004;44(suppl):SP285

Molecular genetic blood group typing by the use of PCR-SSP technique M Prager, BAG Health Care, Lich, Germany Presented at the FDA Workshop on Molecular Methods in Immunohematology Bethesda 25 – 26 Sept 2006 Transfusion 2007;47:54S-59S.

ABO genotyping for diagnosis of unusual ABO blood groups: A Comparative Study in German Blood Donor Centers M Prager, EA Scharberg, FF Wagner, J Burkhart, <u>A Seltsam</u> Presentation at AABB Annual Meeting 2007, October 20 – 23, Anaheim, California, USA



# **BAGenotype Novel Evaluation Software**

K BAGenotype	
Ergebnis       1       2       3       4       5       6       7       8       Genotype       Phanotyp         Monosult       -       +       +       +       +       +       -       A <sup>2</sup> A <sup>2</sup> A <sub>2</sub>	
Probenname: / Name: Sample, S. Proben-ID: / Sample/D: 0000-2008 Geb-Datum: / Birthdate: 01.01.1950 Datum: / Date: 05.05.2008 Unterschrift: / Signature:	Inuckvorschau     X       Inucken     Image: Section of the secti
Kit       00nor/Patient         LOT       0710 AB         Expiry Date       2009-04         Sample.D       0000-2008         Birthdate       01.01.1950         Date       05.05.2008	<text><text><text><text><text><text><text></text></text></text></text></text></text></text>
Fast and easy tool for blood group genotyping using	Elev/CNBAGenetere V1.20 en/CD/inst-hm 20.08.2008

🌯 Start 🚱 🥭 💽 😿 📧 💽

» 💿 Posteinaa..

**BAGene PCR-SSP test kits** 



Bagenot... ■ « ● 18:02

BAGenoty...

🗀 inst

BAGene