

# HISTO TRAY ABC 120/2

Lot. No.: 007T0380

Pos. Well	Serum Lot No.	Anti-HLA	Reaktion Reaction	Pos. Well	Serum Lot No.	Anti-HLA	Reaktion Reaction
HISTO TRAY 120 / 2	1A	105 H 0150	<i>pos. Control</i>	6F	107 H 0120	B 35+53 (51)	
	1B	106 H 0380	<i>neg. Control</i>	6E	105 P 1060	B 35+53+5	
	1C	105 P 0290	B 62	6D	105 H 0810	B 35+53+5+78	
	1D	105 P 0940	B 62+75	6C	103 P 1890	B 37	
	1E	103 P 0480	B 62+63+75+77 (57)	6B	104 P 0280	B 37	
	1F	104 H 1730	B 15+46 (57, 70)	6A	107 H 0320 *	B 60	
	2F	106 P 0730	B 38 (49)	7A	94 P 352	B 60+48	
	2E	105 H 1010	B 16[38+39]	7B	204 H 1250	B 40[60+61]+48 (7)	
	2D	106 H 0200	B 16[38+39]+67	7C	105 P 0600	B 40[60+61]+47+13	
	2C	105 H 1140	B 39+67	7D	103 H 1090	B 41(60)	
	2B	103 H 1170	B 57	7E	106 H 0370	B 41+42+63 (57)	
	2A	106 H 0830	B 17[57+58]	7F	104 H 1560	B 42 (41)	
	3A	107 H 0020	B 17[57+58]+63	8F	204 H 0470	B 47+37+52 (18, 39, 41)	
	3B	106 H 0450	B 18	8E	94 H 402	B 70+21+15+56 (35, 52)	
	3C	103 H 1710	B 18	8D	94 H 400	B 70[71+72]+56 (15,35,21) ↓	
	3D	103 H 1070	B 18 (51)	8C	94 H 403	B70[71+72]+62+35+50+45+41(60)	
	3E	205 H 0090	B 49	8B	105 H 1500	B 73	
	3F	107 H 0070	B49+52	8A	106 H 0310	Cw 1	
	4F	94 P 438	B 21[49+50]+56 (35)	9A	106 H 0290	Cw 2	
	4E	106 H 0690	B 21[49+50]	9B	107 H 0700 *	Cw 3	
	4D	106 H 0160	B 55	9C	105 H 0230	Cw 4	
	4C	107 H 0090	B 55+67	9D	204 H 1230	Cw 5	
	4B	106 H 0570	B 22[54+55+56]+42+7+67	9E	94 H 427	Cw 6	
	4A	106 P 0150	B 56	9F	94 H 428	Cw 6 (4)	
	5A	105 H 0910	B 27+47	10F	94 H 333	Cw 7	
	5B	105 H 1470	B 27	10E	103 H 1410	Cw 8+5	
	5C	106 H 0230	B 27	10D	105 H 0510	Bw 4+A 32	
	5D	105 H 0750	B 27	10C	105 H 0170	Bw 4+B45	
	5E	106 H 0420	B 35	10B	106 H 0190	Bw 6	
	5F	94 P 211	B 35+53 (Cw 4)	10A	104 H 0100	Bw 6+B13	

ABC120/SD/1000/082007

**Bemerkungen/Remarks:**

- wk = schwach/weak
- [ ] = Split-Reaktion/Split-reactions
- \* = Neues Antiserum/new antiserum
- Neue Position/new position
- + e = extra Reaktion/extra reaction
- ( ) = kann Reaktion zeigen/may show reactions
- s.t.wk = manchmal schwach/sometimes weak
- + = monoclonal
- ↓ = carry over possible

## Bw4 und Bw6 Assoziationen Bw4 and Bw6 associations

**Bw4:** B5 · B5102 · B5103 · B13 · B17 · B27 · B37 · B38 (16) · B44 (12) · B47 · B49 (21)  
B51 (5) · B52 (5) · B53 · B57 (17) · B58 (17) · B59 · B63 (15) · B77 (15)

möglich  
possible A9 · A23 (9) · A24 (9) · A2403 · A25 (10) · A32 (19)

**Bw6:** B7 · B703 · B8 · B14 · B18 · B22 · B2708 · B35 · B39 (16) · B3901 · B40 · B4005 · B41 · B42  
B45 (12) · B46 · B48 · B50 (21) · B54 (22) · B55 (22) · B56 (22) · B60 (40)  
B61 (40) · B62 (15) · B64 (14) · B65 (14) · B67 · B70 · B71 (70) · B72 (70) · B73 · B75 (15)  
B76 (15) · B78 · B81

# HISTO TRAY ABC 120/1

## Ergebnisprotokoll / Worksheet

Mikrotestkammern (60) mit vorgetropften HLA-ABC-Antiseren und Kontrollen (human oder monoklonal) zur Typisierung eines Probanden

*Microtest tray (60) with prealiquoted HLA-ABC anti sera and controls (human or monoclonal) for typing one specimen*

Proben Nr./Sample I.D.: ..... Testdatum/Test Date: .....

Name/Name: ..... Untersucher/Tech: .....

Vorname/First Name: ..... Vitalität/Viability %: .....

Geb.-Datum/Birth date: .....

A ..... A ..... B ..... B .....

Cw ..... Cw ..... Bw ..... Bw .....

 7013

**CE 0123**

 007T0380




 2008-08



 ≤ -20°C

Kaninchenkomplement:  
Rabbit Complement:

 007K0721

 2009-06

	Pos. Well	Serum LOT	Anti-HLA	Reaktion Reaction	Pos. Well	Serum LOT	Anti-HLA	Reaktion Reaction
	<b>HISTO TRAY 120 / 1</b>	1A	105 H 0150	pos. Control		6F	304 P 0060	A 30+31
1B		106 H 0380	neg. Control		6E	105 H 0790	A 30+31+B13	
1C		94 P 271	A 1		6D	204 H 0080	A 32	
1D		106 P 0760	A 1		6C	104 H 0310	A 32+25	
1E		107 H 0130	A 1+36		6B	103 H 1540	A 32 (25)	
1F		105 H 0420	A 2		6A	203 H 1050	A 32+74+1+3+11+36	
2F		106 P 0500	A 2		7A	303 H 1060	A 33 (34)	
2E		105 P 0390	A 2+69		7B	106 H 0410	A 33+31 (29)	
2D		204 H 0250	A 3		7C	94 H 434	A 33+10+28	
2C		106 H 0390	A 3		7D	104 P 1240	B 51	
2B		94 P 350	A 3+11		7E	105 H 0800	B 51	
2A		204 P 1200	A 23		7F	105 H 0550	B 5[51+52]+78	
3A		105 P 1480	A 9 [23+24+2403]		8F	105 H 0570	B 5[51+52]+53	
3B		303 H 1040	A 9 [23+24]+80		8E	105 P 0330	B 7+27+73	
3C		204 H 1180	A 24+2403		8D	105 P 0040	B 7 (42)	
3D		106 H 0300	A 25		8C	103 P 0540	B 7+81	
3E		103 P 1570	A 25		8B	94 P 287	B 8	
3F		105 P 0450	A 10[25+26+66+34]		8A	103 P 1690	B 8+59	
4F		94 P 442	A 26+66		9A	106 H 0530	B 8+14+39	
4E		105 H 0350	A 26		9B	107 H 0040	B 44	
4D		103 P 1500	A 11		9C	105 H 0120	B 44+47 (49)	
4C		106 P 0040	A 11		9D	104 P 0800	B 12[44+45]	
4B		204 H 0870	A 28[68+69]		9E	204 H 0500	B 45+82	
4A		303 H 1740	A 28[68+69]		9F	105 P 0010	B 13	
5A		204 H 0730	A 68+34+66 (25, 33)		10F	105 H 0030	B 13	
5B		105 H 0370	A 29		10E	105 H 0020	B 13	
5C		94 H 453	A 29+43		10D	203 P 1730	B 64+8	
5D		103 P 0060	A 29		10C	106 H 0610	B 14[64+65]	
5E		104 H 0750	A 30		10B	105 P 0830	B 14[64+65]	
5F		204 P 0010	A 30		10A	105 H 0900	B 14[64+65]	